

Venue Go Ultrasound System Datasheet

Product Description

Venue Go™ is a take-anywhere ultrasound system that provides the latest technologies to help deliver a simple, fast, and precise solution to the Point of Care ultrasound community. Its innovative design may also be useful for traditional users in a variety of clinical settings.

Venue Go combines a single-surface touch screen interface, which is great for quick clean-up, with cords-off-the-floor docking for transducer protection. It has been designed and tested to allow walk-up simplicity to minimize training. This versatile system allows the user to choose to scan with the Venue Go as a compact system, as a console, or attached to a boom, freeing up much-needed space in an otherwise overcrowded environment. The Venue Go system is built on GE Ultrasound's advanced cSound™ Imageformer providing ploppable image quality for users who want to scan and go. Automated tools help reduce the fatigue of repeated and difficult tasks and help simplify and accelerate clinical workflow. The tap, pinch, and slide gestures used in the Venue Go user interface have the familiar feel of contemporary consumer devices. Underneath it all, Venue Go displays an impressive set of features designed for the novice as well as the demanding user, both of whom need to be ready to scan at any moment.



General Specification

System Dimensions

Height: 348 mm (13.7 in.)

Width: 406 mm (16.0 in.)

Depth: 59 mm (2.3 in.), 127 mm (5.0 in.) with probe holder

Weight: 6.27 kg (13.9 lbs) without probes

System Dimensions with Optional Cart

Height: 1190-1580 mm (46.9 - 62.2 in.)

Width: 505 mm (19.9 in.)

Depth: 480 mm (18.9 in.)

Weight: 32 kg (70 lbs) without probes

Optional Cart Dimensions

Height: 1000-1390 mm (39.4 - 54.7 in.)

Width: 505 mm (19.9 in.)

Depth: 480 mm (18.9 in.)

Weight: 25 kg (55.6 lbs.) without scanner & accessories

Cart design

Quick release / mount of scanner to cart

Multipurpose handle - customizable

Probe-holder insert

Gel cup

Barcode scanner holder

Built-in AC adapter

Horizontal swivel of display on cart: 230° (115° to each side)

Backward tilt of display on cart: 40°

Forward tilt of display on cart: 10°

Up / down adjustment range of display on cart: 390 mm (15.5 in.)

Vesa mount: 75x75 mm (2.98 x 2.98 in.)

Four 5" swivel-casters with brake & swivel lock on each caster

Cradle may be provided as an accessory

Anti-theft cable attachment area for optional cable

Electrical Power

Voltage: 100-240 V AC

Frequency: 50/60 Hz

Power: Max. 250 VA

Standard lithium ion battery with up to two hours scan time when fully charged

Recharge time: 2 hours

No. of charge cycles to age 70% capacity = 300

Battery power breaker switch for transportation

Design

Compact

Full touch user interface

3 active probe ports

Multipurpose handle, for probes, gel, handle or accessories.

cSound™ beamforming technology

CPU – Skylake-Intel® 6th Generation, i7 quad core

Operating System – Microsoft® Windows® 10

Single Integrated Solid State Drive (128 GB)

Archive size: > 60 GB

AC adapter

Adjustable rear support stand for use on flat surfaces

System Connectivity

USB 3.0 4 connectors

Ethernet LAN 10/100/1000 Base T

HDMI Output available for compatible devices

Wireless LAN capability, with WiFi strength indication symbol

User Interface

Touch Screen

Multi-touch user-interface with gesture recognition

Probe & application quick select

Home screen for easy workflow

Clinical applications oriented workflow

Mode specific controls

Alphanumeric touch-keyboard in selected language

Touch Screen (continued)

Measurements (full touch)

Measurement worksheet

Thumbnails and image viewer

Annotations and findings

Body marks

Utility settings

Patient information

User interface in 13 configurable languages

Touch gesture recognition

The following controls can be adjusted by intuitive gestures

Gain on all scanning modes

Depth of B, color and M-Mode

Focus position

Scale, baseline and sweep-speed of PW or CW modes

Sweep speed when in M-Mode or Anatomical M-Mode

Zoom in/out

Cine review

Color ROI position/dimension

PW cursor and SV gate length & location

M-Mode and anatomical M-Mode cursor adjustment

Center line on/off

Display

15.6 in. Multi-touch High-Resolution Color LCD, 16:9 ratio

Resolution: 1920 x 1080 pixels

Horizontal/Vertical viewing angle: +/-80°

Display Frame rate (Hz): 60Hz

Light sensor for Auto ambient light compensation

Integrated speaker

Audio Volume Adjustment

Battery and wall/AC power indicator displays even when system is turned off

Battery time remaining indicator

3 probe holders and multipurpose holder on display

Flexible display positions

Display (continued)

Screen cleaning mode capability

Sealed display

Diagnostic LCD

Capacitive sealed power button

LED Color bar for system identification and Low Battery warning

System Overview

Applications

Abdominal

Cardiac (Adult and Pediatric)

Thoracic/Pleural

Neonatal (head / abdomen / hip)

Adult and neonatal cephalic

Breast

Gynecology

Musculoskeletal

Obstetrical

Nerve

Interventional Procedure

Pediatric

Small Parts

Urological

Vascular

Vascular access

Transvaginal

Transrectal

Intraoperative

Ophthalmic (optional)

Cardiac TEE

Transducer Types

Convex Array: C1-5-RS

Microconvex Array: 8C-RS, E8C-RS

Phased Array: 3Sc-RS, 6S-RS

Linear Array: 9L-RS, 12L-RS, L4-12t-RS, L8-18i-RS

TEE Adult Probe: 6Tc-RS

Operating Modes

B-Mode

Coded Harmonic Imaging (CHI)

M-Mode

Anatomical M-Mode

Color Flow Mode (CFM)

M-Color Flow Mode

Power Doppler Imaging (PDI)

Continuous Wave Doppler (CWD)

Pulse Wave Doppler (PWD) with high PRF

Tissue Doppler Imaging (TDI)

Standard Features

DICOM® 3.0 Connectivity

Integrated Solid State Drive

Automatic Optimization

CrossXBeam™ Compounding

Speckle Reduction Imaging (SRI-HD)

Virtual Convex

Virtual Apex

Raw Data Processing

Real-time automatic Doppler calculations

On-board User Manual

Electronic Documentation eIFU access through the internet

Cine loop storage and playback

InSite™ ExC capability

Patient Information Entry

Measurement & Annotation Package

Editable calculations

Cardiac Calcs

OB Calcs and Tables

Multi Gestational Calcs

Fetal growth trending

Musculoskeletal Calcs

Gynecological Calcs

Urological Calcs

Full Screen mode

Standard Features (continued)

Split Screen

Needle Recognition

Zoom

Patient data encryption

Patient data anonymization

User-defined Presets

Lung diagram

Lung scoring

Software Options

WiFi Activation

Ophthalmic

Hardware Options

USB thermal B&W printer

USB ECG with 3 electrodes (AHA/IEC), including External ECG interface

ECG module is type CF Defibrillator-Proof

USB Memory Stick

Internal WiFi (Optional) supports WLAN, WPA/WPA2 LEAP/PEAP, FIPS-2 WiFi, WPA- Enterprise and WPA2 Enterprise protocols

WiFi bandwidth is 2.4GHz to 2.4835GHz and 5.15GHz to 5.85GHz, dependent on country.

Bar code reader with multi-field mapping support

Large, removable storage bin (Standard on cart)

Storage bin dimensions: L: 37cm, W 15cm, H: 18cm, (14.7 x 6.0 x 7.1 in.)

Mountable intracavitary probe holder

Disposable sterile touch-screen cover

Anti-theft cable

Advanced Applications

Venue Go was designed to host a variety of advanced clinical applications dedicated to specific point-of-care uses.

To do this, we use Machine Learning (ML) and Artificial Intelligence (AI) to develop algorithms that are then used in the system to enable the shift of cognitive burden from ultrasound acquisition to clinical reasoning at the point of care.

AI Enabled Auto Tools

Auto VTI: Calculates the velocity time integral (VTI), stroke volume, CO Flux and cardiac output in a single step. Like the other tools, it includes a quality indicator to assist with image acquisition.

ML: Finding the right location and tracing the spectral doppler signal and calculate and display the VTI while you scan.

AI: AI helps localize the acquisition.

Auto B-Lines: Highlights and counts B-lines in real-time. Hit freeze and Venue Go displays the frame with the highest B-line count.

ML: Automatically distinguishes between real B-lines and all other artifacts is hard to "program", but ML can enable it.

AI: Classification between artifacts vs B-lines.

Auto-IVC: Measures IVC collapsibility. IVC diameter changes (Collapsibility or distensibility index) are measured and displayed in real-time upon completion of each respiratory cycle.

ML: Auto IVC: One respiratory cycle Distensibility Index (DI) or Collapsibility Index (CI) can be displayed, and calculated.

AI: Identify the location and position IVC for measurements.

Lung Tool: See all ultrasound lung findings in one view. Keeps track of your segmental lung assessment. Like the other tools, it can be helpful in showing trends in response to therapy.

VTI Trending: Allows you to quickly visualize the trend and help you determine a next course of action in treatment.

eFAST diagram: Allows scanning, one-tap allocation and quick review of images and findings belonging to different zones of the eFAST and FAST exam.

Display Modes

Simultaneous Capability (duplex or triplex), probe dependent

- B/PW
- B/CW
- B/TDI
- B/CFM
- B/PDI
- B/M
- B/AMM
- Dual B + CFM
- Dual B + PDI
- Needle Recognition + CFM or PDI
- B/CFM/PW or B/CFM/CW

Selectable Alternating Modes, probe dependent

- B/M
- B/PW
- B/CW
- B/TDI
- B/AMM
- B/M/CFM
- B/M/PDI
- B/CFM/PW
- B/CFM/CW
- B/PW/PDI

Multi-image Split Screen

- Live and/or frozen
- Independent Cine playback
- ECG-synchronized Cine playback

Zoom

- Pan zoom (read zoom)
- Expand the image to full screen

Colorized Image

- Colorized B
- Colorized B+M
- Colorized PW
- Colorized CW

Timeline Display

Independent Dual B/PW/CW Display

Display formats: Top/bottom or side-by-side

Format Size: Vert1/3 B, Vert1/2 B, Vert2/3 B, Horiz 1/2B, TL Only

Format, switchable after freeze

Security

Encryption for internal data

Encryption for external media

DICOM® TLS

Whitelisting

LDAP interface for user credentials

Customized policies for user accounts

Security audit trail - for relevant security events

Display Annotation

Institution/Hospital Name

Date: 2 types selectable MM/DD/YY, DD/MM/YY

Time: 2 types selectable 24 hours, 12 hours

Operator Identification

Patient Name: Last, First

Patient ID

Imaging parameters

Gestational Age from LMP/EDD/GA

Patient Age

Power Output Readout:

- MI: Mechanical Index
- TIS: Thermal Index Soft Tissue
- TIC: Thermal Index Cranial (Bone)
- TIB: Thermal Index Bone

TEE probe: Angle of rotation and temperature

System Status (real-time, cine-played or frozen)

Probe Orientation Marker: Coincides with orientation marking on the image monitor

Measurement Results Window: presettable display location

Probe Type and symbol

Application Name

Annotation Library: up to 21 preset labels (depending on application)

Customizable annotations: more than 24 available for each application

Comments available in Live scan mode and Freeze mode

Body marks available for each application

Arrows available in Live scan mode and Freeze mode

Center-line marker on/off on convex and linear probes

Heart rate – real time from ECG Signal

Battery status and remaining live scanning time indicator

Anatomy specific presets

Image Preview

Gray/Color Bar

Measurement Summary Window

Imaging Parameters by Mode (current mode/see below)

Focal Zone Markers

Image Management Menu: Menu, Delete, and Image Manager

Image Palette

System Parameters

System Setup

User Programmable Preset Capability

Factory Default Preset Data

Factory Default Application Data

User Interface Languages: Brazilian Portuguese, Danish, Dutch, English, Finnish, French, German, Italian, Norwegian, Portuguese, Russian, Spanish and Swedish

User Manual Languages: Brazilian Portuguese, Bulgarian, Chinese, Croatian, Czech, Danish, Dutch, English, Estonian, Finnish, French, German, Greek, Hungarian, Indonesian, Italian, Japanese, Kazakh, Korean, Latvian, Lithuanian, Norwegian, Polish, Portuguese, Romanian, Russian, Serbian, Slovakian, Slovenian, Spanish, Swedish, Turkish, Ukrainian and Vietnamese

Operation Error Message Display

System boot from full shutdown: < 75 sec.

Sleep / Standby fast boot up: < 7 sec. on average

Probe Loading: < 3 sec.

Setup wizard for easy set up

Image Processing and Parameters

Digital beamformer based on cSound™ Architecture

Displayed Imaging Depth

- Minimum Depth Range of Field: 0 - 1.0 cm (probe dependent)
- Maximum Depth Range of Field: 0 - 36 cm (probe dependent)

Continuous Dynamic Receive Focus/Aperture

Multi-Frequency/Wideband Technology

Reporting

Qview is an application running on Venue Go that securely provides web access to Telexy Healthcare's Qpath product and enables exam reporting at the bedside. The exam images and report worksheet are displayed on the Venue Go monitor just as they would be if you logged into Qpath from a computer workstation. Qpath is not available in all countries, please verify with your GE representative. Qpath is an application under license from Telexy Healthcare.

Cine Memory

Able to save an image or cine loop during scanning

Configure retrospective or prospective cine loop storage

Configurable cine loop length (up to 180 sec. or 100 heart-cycles)

Image Archive and Connectivity

Image Browser: Previewing of previous archived images as well as current stored patient images

View images side-by-side from different exams

Image Management (USB memory-stick or removable hard-disk drive)

Delete selected exams or images

Review in Full Image Area

Direct print to approved printer

Live Scan "Store": Configure Store button to save an image during live scanning

Archiving Format

JPEG

MPEG4 /H.264

AVI, DICOM, Raw Data+DICOM, WMV

Capture Area

Image Area

Full Screen

DICOM

DICOM store

Worklist query

Multi-frame DICOM

Network store to a shared folder

Scanning Parameters

B-Mode

Brightness mode. Real time displays of a two dimensional cross section of a three-dimensional soft tissue structure. Ultrasound echoes of different intensities are mapped to different gray scale or color values in the display.

B/M Acoustic Output: 0–100%

Virtual convex: On/Off on Linear array probes

Image Reverse: On/Off

B Colorize maps: 4 types

B-Mode (continued)

Thermal Index: TIC, TIS, TIB

CrossXBeam Compounding

TGC plot and 8 graphical TGC sliders

"Auto" image optimize

Width (FOV)

Speckle Reduction SRI HD: 5 increments

Edge Enhance: 5 increments

Probe dependent focal zones

Line Density: up to 4 increments, application dependent

Frame Average: up to 8 increments, application dependent

Angle (degree): application dependent, 10°– 133°

Gray Scale Map: up to 11 types, application dependent

Gain and Depth Adjustment

CHI (Coded Harmonics Imaging)

Dynamic Range: 36–102 dB, application dependent

Depth range: 0 - 36 cm, application dependent

Rejection: 6 increments

Color Flow Mode (CFM) or Color Doppler

A real-time two-dimensional cross-section image of blood flow.

Color gradient used to represent directional blood flow (velocity, variance, power and/or direction) prioritized over amplitude.

Steer

Compression

Base Line: application dependent

Invert: On/Off

Focus Depth

Power - Acoustic Output: 0–100 %

Line Density: 3 increments

Frame Average: up to 11 increments, application dependent

PRF: 0.3K–22 KHz: probe dependent

Sample Volume

Gain: 0–40 dB

Wall Filter: probe dependent

Lateral and radial filters

Vertical Size (mm): default presettable

CF/PDI Center Depth (mm): default presettable

Color Flow Mode (CFM) or Color Doppler (continued)

Color Frequency

Focal Number: 1

Color Maps: up to 6 types

Tissue Priority: 13 steps

Power Doppler Imaging (PDI) Mode

Color gradient used to represent blood flow using amplitude vs. velocity shift (Color Doppler). Prioritizes amplitude over direction.

ROI Position

ROI Size

Focus Position

Steer

Compression: 15 increments

Line density: 3 increments

PDI Maps: 5 types

Acoustic Output: 0–100%

Packet Size: 8–24, probe dependent

Spatial Filter: 6 steps

Frame Average: 7 steps, probe dependent

PRF: 0.3K–11.4K Hz: depth dependent

Power Threshold: 10–100%

CF/PDI Vertical Size: default presettable

CF/PDI Center Depth: default presettable

CF/PDI Focal Number: 1

Gain: 0–40 dB

Wall Filter: 4 increments, probe dependent

CF/PDI Frequency: 2–4 increments, probe dependent

PDI Flash suppression: 7 increments

PDI Radial & Lateral averaging: 5 increments

PDI Tissue priority: 10 increments

M-Mode/Anatomical M-Mode

Motion mode. Soft tissue structure is presented as scrolling display, with depth on the Y-axis and time on the X-axis. Anatomical M-Mode (AMM) Allows M-Mode on stored 2D cine clip. Facilitates arrhythmia assessment and cardiac measurements.

Sweep Speed

M Color

MM Display Format: Vert1/3 B, Vert1/2 B, Vert2/3 B, Horiz 1/2B, TL Only

B/M Acoustic Output: 0–100%

Rejection: 6 increments

Dynamic Range: 36–96 dB

Edge Enhance

Gray Scale Maps

M Gain: +/- 20 dB delta from B

PW/CW Mode

Pulse Wave Doppler (PW), Continuous Wave Doppler (CW) are used for displaying the speed, direction, and spectral content of blood flow at selected anatomical sites.

PW velocity

- Max: 870 cm/s (19,800 Hz)
- Min: 3 cm/s (140 Hz)

CW velocity

- Max: 1,310 cm/s (40,000 Hz)
- Min: 5cm/s (260 Hz)

Gray Scale Maps

Compression: 12 steps

Base Line: 5–95%

SV Gate size: 1 - 16 mm, application dependent

Angle Correction: +/- 90°

Spectral Color: 6 types

PW Sweep Speed

PW/CW Display Format: Vert1/3 B, Vert1/2 B, Vert2/3 B, Horiz 1/2B, TL Only

PW Acoustic Output: 0–100%

Adjustable cursor steering and gate position

Selectable real-time duplex and triplex operation

Spectral Averaging

PW/CF Ratio: 1, 2, 4

PW/CW Mode (continued)

Rejection

Doppler Gain: 0–32 dB

Wall Filter: probe dependent

PW Angle Steer: 0 +/- 20°

Sample Volume Depth: 29 increments default presettable

Audio Volume

PW Frequency 2–4 steps: probe dependent

M-Color Flow Mode

Overlays color on the M-Mode trace

Coded Harmonic Imaging (Tissue Harmonics) (CHI)

Enhances near field resolution for improved small parts imaging as well as far field penetration. Diminishes low frequency amplitude noise and provides clarity to needle, anatomy and motion.

TVD (TDI)

Tissue Velocity Doppler provides tissue spectral information for selected Doppler sample.

Needle Recognition Mode

Provides accurate display of the needle, anatomy and motion even in Color and Power Doppler.

Includes ability to adjust needle gain and angle.

Available on all linear and some convex probes.

Automatic Optimization

Continuous Tissue Optimization

Auto Spectrum Optimization

CrossXBeam Compounding

Provides 3 angles of spatial compounding

Compatible with Color Mode

SRI-HD

Speckle Reduction Imaging provides multiple levels of speckle reduction

Compatible with/ B-Mode, Color

Virtual Convex

Provides a trapezoidal field of view on all linear array probes, for wider FOV in the far field

Virtual Apex

Provides wider FOV in the near field

Supported on some presets on Sector probes

Measurements and Calculations

B-Mode Calcs

Distance

Circumference (Ellipse/Trace)

Area (Ellipse/Trace)

% Stenosis

Angle between 2 lines

Ratios

Depth from Probe Surface

M-Mode Calcs

Distance

Time

Slope

Heart Rate

Doppler Calcs

Velocity

Frequency

Time

Acceleration

Heart Rate

Auto Doppler Trace function with automatic calcs

Time averaged max/mean velocity

Ratios

PI (Pulsatility Index)

RI (Resistivity Index)

Vascular Measurements and Calcs

Upper/Lower

Artery/Vein

Summary Worksheet

Obstetrics Measurements and Calcs

Gestational Age Calculation

Multi-Gestational Calculation, up to 3 fetuses

EFW Calculation

Summary Worksheet

Fetal Trend Graph

Abdominal Circumference (AC)

Amniotic Fluid Index (AFI)

Area

Diameter (APTD-TTD)

Bi-parietal Diameter (BPD)

Crown Rump Length (CRL)

Estimated Fetal Weight (EFW)

Femur Length (FL)

Gestational Sac (GS)

Head Circumference (HC)

Humerus Length (HL)

Occipito frontal Diameter (OFD)

Cardio-Thoracic Area Ratio (CTAR)

Fetal Trunk Cross-Sectional Area (FTA)

Spine Length (SL)

Comparison of multiple fetus data on a graph and a worksheet

Gynecology Measurements and Calcs

Ovarian Follicle Measurements

Summary Worksheet

Urology Measurements and Calcs

Volume Measurements

Summary Worksheet

Musculoskeletal Measurements and Calcs

Labeled measurements

Cardiac Measurements and Calcs

Ventricle, Atrium, Valve Measurements

Summary Worksheet

Probes (all optional)

3Sc-RS Wide Band Phased Array

Applications: Abdominal, Adult Cardiac, Adult Cephalic, Gynecology, Obstetrics, Ophthalmic, Pediatric Cardiac, Pediatric/Neonatal, Peripheral Vascular, Thoracic/Pleural, Urology, Vascular Access

Imaging Frequency B-Mode

Fundamental: 2.5 - 4.5 MHz

Harmonics: 1.6/3.2 - 1.9/3.7 MHz

CFM Imaging Frequency: 1.9 - 3.6MHz

6S-RS Wide Band Phased Array

Applications: Abdominal, Adult Cardiac, Gynecology, Neonatal Cephalic, Pediatric Cardiac, Pediatric/Neonatal, Peripheral Vascular, Thoracic/Pleural, Urology, Vascular Access

Imaging Frequency B-Mode

Fundamental: 4.0-6.5 MHz

Harmonics: 2.7/5.4 - 3.2/6.4 MHz

CFM Imaging Frequency 2.6-5.9 MHz

8C-RS Wide Band Microconvex

Applications: Abdomen, Adult Cardiac, Gynecology, Neonatal Cephalic, Nerve Block, Pediatric Cardiac, Pediatric/Neonatal, Peripheral Vascular, Thoracic/Pleural, Urology

Imaging Frequency B-Mode

Fundamental: 6.0 - 10.0 MHz

Harmonics: 3.5/7.0 - 5.0/10.0 MHz

CFM Imaging Frequency: 4.3 - 5.6 MHz

C1-5-RS Wide Band Convex

Applications: Abdominal, Adult Cardiac, Fetal/Obstetrics Gynecology, Musculoskeletal Conventional, Musculoskeletal Superficial, Nerve Block, Pediatric/Neonatal, Pediatric Cardiac, Peripheral Vascular, Thoracic/Pleural, Urology, Vascular Access

Imaging Frequency B-Mode

Fundamental: 3.0 - 5.0 MHz

Harmonics: 1.5/3.0 - 2.5/5.0 MHz

CFM Imaging Frequency: 1.8 - 3.6 MHz

E8C-RS Wide Band Microconvex Intracavitary

Applications: Gynecology, Obstetrics, Transrectal, Transvaginal, Urology

Imaging Frequency B-Mode

Fundamental: 6.0 - 10.0 MHz

Harmonics: 3.5/7.0 - 5.0/10.0 MHz

CFM Imaging Frequency: 4.3 - 5.6 MHz

L4-12t-RS Wide Band Linear with Buttons

Applications: Abdominal, Adult Cardiac, Gynecology, Musculoskeletal Conventional, Musculoskeletal Superficial, Nerve Block, Ophthalmic, Pediatric/Neonatal, Pediatric Cardiac, Peripheral Vascular, Small Parts, Thoracic/Pleural, Urology, Vascular Access

FOV (max): 38.4mm

Imaging Frequency B-Mode

Fundamental: 7.0 - 12.0 MHz

Harmonics: 3.5/7.0 MHz - 5.5/11.0 MHz

CFM Imaging Frequency: 4.5 - 8.3 MHz

4 Configurable Buttons to Support Various Operations

Biopsy Available: Civco Infiniti Plus™ style bracket. (P/N: H48392LT / 742-429).

Civco AccuSITE™ Out-of-Plane Ultrasound Needle Guide (P/N): H48392LL / 742-359-10)

12L-RS Wide Band Linear Array

Applications: Abdominal, Adult Cardiac, Gynecology, Musculoskeletal Conventional, Musculoskeletal Superficial, Neonatal Cephalic, Nerve Block, Ophthalmic, Pediatric/Neonatal, Pediatric Cardiac, Peripheral Vascular, Small Parts, Thoracic/Pleural, Urology, Vascular Access

FOV (max): 38.4 mm

Imaging Frequency B-Mode

Fundamental: 7.0 - 12.0 MHz

Harmonics: 3.5/7.0 MHz - 5.5/11.0 MHz

CFM Imaging Frequency: 4.5 - 8.3 MHz

Biopsy Available: Civco Infiniti Plus™ style bracket (P/N: H48392LT / 742-429).

Civco AccuSITE™ Out-of-Plane Ultrasound Needle Guide (P/N): H48392LL / 742-359-10)

L8-18i-RS Wide Band Linear hockey-stick shape

Applications: Abdominal, Gynecology, Intraoperative, Musculoskeletal Conventional, Musculoskeletal Superficial, Neonatal Cephalic, Nerve Block, Pediatric/Neonatal, Peripheral Vascular, Small Parts, Thoracic/Pleural, Urology, Vascular Access

FOV (max): 25.2 mm

Imaging Frequency B-Mode

Fundamental: 9.0 - 16.0 MHz

Harmonics: 5.0/10.0 - 9.0/18.0 MHz

CFM Imaging Frequency: 5.6-12.5 MHz

9L-RS Wide Band Linear Array

Applications: Abdominal, Adult Cardiac, Fetal/Obstetrics, Gynecology, Musculoskeletal Conventional, Musculoskeletal Superficial, Neonatal Cephalic, Nerve Block, Ophthalmic, Pediatric/Neonatal, Pediatric Cardiac, Peripheral Vascular, Small Parts, Thoracic/Pleural, Urology, Vascular Access

FOV (max): 44.16 mm

Imaging Frequency B-Mode

Fundamental: 3.0 - 8.0 MHz

Harmonics: 2.8/5.6 MHz - 4.5/9.0 MHz

CFM Imaging Frequency: 3.1 - 6.3 MHz

Biopsy Available: Civco Infiniti Plus™ style bracket (P/N: H45281BL / 742-416)

6Tc-RS Wide Band Phased Array Transesophageal (TEE)

Applications: Cardiac-TEE

Fundamental: 3-8 MHz

Harmonics: 3.0/6.0 MHz

CFM Imaging Frequency: 3.1-6.3 MHz

Note About Sales in your Country

Not all features, products, probes, or peripherals described in this data-sheet may be available or cleared for sale in all markets. Please contact your local GE Ultrasound representative to get the latest information.

Safety Conformance

Approved for ANSI/AAMI ES60601-1 by Nationally Recognized Testing Laboratories

Approved for CAN/CSA-C 22.2 No.60601.1 by Nationally Recognized Testing Laboratories

CE Marked to Council Directive 93/42/EEC on Medical Devices

Compliant with DIRECTIVE 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) requirement.

- IEC¹ 60601-1 Safety requirements for medical electrical systems
- EC¹ 60601-1-2 General requirements for safety – Electromagnetic compatibility
- IEC¹ 60601-1-6 Medical Electrical Equipment— General requirements for basic safety and essential performance— Usability
- IEC¹ 60601-2-37: Particular requirements for the safety of ultrasonic medical diagnostic and monitoring equipment
- ISO 10993-1: Biological evaluation of medical devices (applicable to the probes)
- NEMA UD3 Acoustic output display (MI, TIS, TIB, TIC)

¹Including national deviations

About GE Healthcare

GE Healthcare provides transformational medical technologies and services to meet the demand for increased access, enhanced quality and more affordable healthcare around the world. GE (NYSE: GE) works on things that matter - great people and technologies taking on tough challenges. From medical imaging, software & IT, patient monitoring and diagnostics to drug discovery, biopharmaceutical manufacturing technologies and performance improvement solutions, GE Healthcare helps medical professionals deliver great healthcare to their patients.

GE Healthcare
9900 Innovation Drive
Wauwatosa, WI 53226
U.S.A.

www.gehealthcare.com

Imagination at work

© 2019 General Electric Company – All rights reserved.

Product may not be available in all countries and regions. Full product technical specification is available upon request. Contact a GE Healthcare Representative for more information. Please visit www.gehealthcare.com/promotional-locations.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, GE Monogram, cSound, CrossXBeam, InSite, and Venue Go are trademarks of General Electric Company or one of its subsidiaries. GE Healthcare, a division of General Electric Company. DICOM is the registered trademark of the National Electrical Manufacturers Association. Intel is a trademark of Intel Corporation or its subsidiaries in the U.S. and/or other countries. Infiniti Plus and AccuSITE are trademarks of Civco Medical Solutions. Microsoft and Windows are trademarks of Microsoft Corporation. All other third party trademarks are the property of their respective owner. GE Medical Systems, Inc., doing business as GE Healthcare.



DOC2090867
June 2019